
## INVACARE® MANUAL BED





The **Invacare® Manual Bed** offers the same convenient features as electric models, but adjusts manually for greatest economy. Manual cranks located at the foot end of the bed raise and lower upper body and knees, as well as bed height.

Invacare beds, mattresses and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit **www.invacare.com/bedsafety** for additional safety and entrapment information.

## **Features**

Patented universal bed ends simplify inventory – store just one universal bed end instead of two different bed ends

Patented universal bed ends can be used as either the headboard or footboard on the following Invacare bed models: 5307IVC, 5310IVC, 5410IVC, 5410LOW and VC5310

High-impact bed end panels are stronger than fiberboard and have better impact and scratch resistance

Two locking and two non-locking casters with non-marring rubber tread for secure positioning

Enhanced caster sockets provide better retention during delivery and setup

## **HCPCS Code: E0293**





1. Patented universal bed ends















Model #	Dimensions	Sleep Surface	Weight Capacity	Shipping Weight	Carton Dimensions	Regulation Compliance	Limited Warranty
5307IVC	88" L 15"-23" H 36" W	36" W	Safe Working Load: 450 lb.* 1ax. Patient Weight 350 lb.		N/A	UL962 CSA Certified	5 Years Welds d 2 Years Mechanical
5000IVC				45 lb.	43" L		5 Years Welds
Headspring					8" W		2 Years
					36" D		Mechanical
5770IVC				62 lb.	47" L		
Footspring					9" W		
					36" D		
5301IVC	•	•		55 lb.	39" L		
Universal Bed End					7" W		
					24" D		

<sup>\*</sup>Includes patient, mattress and all accessories.



Invacare Corporation www.invacare.com

USA One Invacare Way Elyria, Ohio 44035-4190

Unit 8

©2015

(800) 333-6900 **Canada** 

570 Matheson Blvd. E.,

Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Specifications are subject to change without notification.
Form No. 11-215
Rev. 08/15
1506103

Yes, you can.